

3.0mm

Anchor with  
one suture



Anchor with One #2 Suture


Tool Used for  
Implantation

**DFX0028A**

2.3mm or 2.5mm drill  
bit with positioning tool



**ANCHOR IMPLANTED  
BY PUSH-IN METHOD**



**SUGGESTED USE**  
Shoulder-Labral Repair

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)



3.5mm

Anchor with  
Two Sutures




Anchor with Two #2 Sutures


Tool Used for  
Implantation

**DFX0029A**

2.8mm or 3.0mm drill  
bit with positioning tool



**ANCHOR IMPLANTED  
BY PUSH-IN METHOD**



**SUGGESTED USE**  
Shoulder-Labral Repair

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)



4.5mm

Knotless Anchor  
with One Suture




Tool Used for  
Implantation

**DFX0030A**

4.5mm Punch



**ANCHOR IMPLANTED  
BY PUSH-IN METHOD**



**SUGGESTED USE**

Shoulder - Lateral Row, RCR

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)



4.75mm

Double-loaded  
Threaded Anchor



Anchor with Two Fiber-Connect Tapes

Tool Used for  
Implantation


**DFX0043A**

Universal Punch  
DM-4818

4.5/4.75mm Punch Tap



**ANCHOR IMPLANTED  
BY THE SCREW-IN METHOD**



**SUGGESTED USE**

Shoulder-Medial Row, RCR

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)



# 5.5mm

## Double-loaded Threaded Anchor



Anchor with Two Fiber-Connect Tapes

Tool Used for  
Implantation


**DFX0040A**

Universal Punch  
DM-4818

5.5mm Punch Tap



**ANCHOR IMPLANTED  
BY THE SCREW-IN METHOD**



### SUGGESTED USE

Shoulder -Medial Row, RCR

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)



# 5.5mm

Triple-loaded  
Threaded Anchor



Anchor with Two Sutures & One Tape

Tool Used for  
Implantation

DFX0036A

Universal Punch  
DM-4818

5.5mm Punch Tap



**ANCHOR IMPLANTED  
BY THE SCREW-IN METHOD**



### SUGGESTED USE

Shoulder-Medial Row, RCR

This anchor is approved by FDA for the following applications:  
Shoulder, Elbow, Knee, Foot and Ankle  
(Refer to IFU for details)

