

# PURmarrow360™

Optimal Recovery of Mesenchymal Stem Cells

## SAFER

Remove concern of disease transmission  
Superior to cell-based matrix

## FASTER

Reduce time and personnel cost  
No centrifugation needed

## BETTER

High cell colonizing units  
Optimize collection of MSCs  
Prevent peripheral blood collection

## REGENERATION

Induce osteogenesis and tissue healing



THE FUTURE OF REGENERATIVE MEDICINE™

[www.PURbiologics.com](http://www.PURbiologics.com)


A HippoFi™ Company

# PURmarrow360™

## Optimal Recovery of Mesenchymal Stem Cells

### EXCLUSIVE DESIGN

- Closed and sealed distal tip
- Patented positive-portal™




### OPTIMAL RECOVERY

- High cellular yield, requiring minimal marrow aspiration
- Safer cell recovery, without risk of cellular sheering


## PUR® vs COMPETITORS

### Initial Volume Aspirated

7-10ml




60ml



### Final Volume

7-10ml



7ml



High cellular yield to volume

Remove concern of disease transmission

Reduce personnel cost

Optimize osteogenesis




### Reduce OR Time

1-2 min.




15-20 min.



PRODUCT NUMBER

PURmarrow-360

ORDERS@PURbiologics.com



PUR®  
biologics

www.PURbiologics.com

A HippoFi™ Company